

tem®

Redefining Bleeding Control



ROTEM® *delta*

ROTEM® *platelet*

Targeted therapy stops the bleeding.

Easy and safe handling
Fast therapeutic decisions
Advanced diagnostic safety



ROTEM® system in clinical use




Ensuring both, the quality and the supply of blood products is becoming increasingly difficult. This is leading to a more rational use of the limited resource “blood” and to a critical view on the need of blood transfusions.

Blood loss, anaemia and blood transfusions are independent predictors for worse outcomes and patients’ quality of life. This includes an increase in morbidity and mortality, as well as a prolongation of the average length of hospital stay. In short, transfusions typically lead to avoidable complications and costs.

In this context, the targeted ROTEM®-based bleeding control solution is integral to Patient Blood Management that is predicated on both pre-emptive and reactive blood-saving measures. In fact, WHA 63.12, all 193 WHO member states have been asked to implement the concept of PBM in a timely manner.

ROTEM®- based bleeding control



ROTEM® *delta* haemostasis analyser



Complicated bleeding situations can occur intra- and post-operatively. They can be life threatening and always require immediate action. A fast differential diagnosis is vital. It is also the basis of a targeted therapy.

The ROTEM®- analysis offers reliable results within 5-10 minutes and provides critical information about the efficacy of the therapy. Additionally the ROTEM®- analysis enables the monitoring and the course-modification of the therapy as needed.

The result is significantly improved patient outcome and lower healthcare cost per episode with subsequent cost savings benefit.

A comprehensive reagent portfolio provides a differential diagnosis

Liquid reagents	in-tem®	r ex-tem®	fib-tem®	ap-tem®/ t ap-tem®	hep-tem®
Single use reagents	in-tem® S	ex-tem® S	fib-tem® S	ap-tem® S	hep-tem® S
	Fast assessment of clot formation, fibrin polymerization and fibrinolysis via the intrinsic pathway	Fast assessment of clot formation, fibrin polymerization and fibrinolysis via the extrinsic pathway	Fast analysis without platelets; qualitative assessment of fibrinogen status	Fast detection of lysis when compared with EXTEM via fibrinolysis inhibition	Specific detection of heparin when compared with INTEM via heparin neutralisation



The ROTEM® *delta* haemostasis analyser measures kinetic changes of the clot elasticity of whole blood samples. It allows quantitative and qualitative assessment by measuring different parameters of the clot status of the blood sample. A comprehensive set of assays permits a differential diagnosis.

The differential ROTEM®- analysis can be performed at the patient's point of care and provides information about hyperfibrinolysis, dilutional coagulopathies, substitution of fibrinogen, factors or platelets as well as the control of heparin or protamine dosage.

ROTEM® *platelet* Proven technologies in one system



The ROTEM® *platelet* measures platelet aggregation in whole blood samples using impedance aggregometry. The new device is run in conjunction with the ROTEM® *delta*, and is compatible with all existing ROTEM® *delta* models with serial numbers > 2000. The ROTEM® *platelet* characteristics are:

- Dedicated single use cuvettes with electrodes
- 2 channels can be used simultaneously
- Measurement time is 6 min, while allowing for emergency samples
- 3 different parameters are available:
 - A6 (amplitude at 6 min in Ohm)
 - MS (maximum slope of the aggregation graph in Ohm/min)
 - AUC (area under the curve in Ohm*min)

Measurements can be performed on the ROTEM® *platelet* device while running measurements on the ROTEM® *delta* system simultaneously.

Single use reagents

adp-tem®

e.g. for the detection of ADP receptor blockage (Clopidogrel)

ara-tem®


e.g. for the detection of cyclooxygenase inhibitors (Aspirin®)

trap-tem®

e.g. for the detection of GP IIb/IIIa receptor antagonists (Abciximab)



ROTEM® *platelet* detection principle




Easy and safe handling



- Single use reagents for fast and reliable results
- Automated pipette for standardised volumes
- Easy operation via touch screen
- Graphical, step-by-step instruction that simplify the test performance
- Integrated learn programme with treatment algorithms and case reports from experts
- Integrated troubleshooting for fast help


Graphical instructions with touch screen



Put the pipette tip into the gently mixed citrated blood sample and confirm

EXTM S	INTEM S	FIBTEM S	APTEM S
1: Huber, Andrea ST: 15:49:21 RT: 00:01:46 CT: 52 s [0038 -- 0079] CFT: --- α: --- A10: ---	1: Huber, Andrea ST: 15:50:07 RT: 00:01:00 CT: * 29 s [0100 -- 0240] CFT: --- α: --- A10: ---	1: Huber, Andrea ST: 15:50:40 RT: 00:00:26 CT: --- CFT: --- α: --- A10: ---	1: Huber, Andrea ST: --- RT: --- CT: --- CFT: --- α: --- A10: ---

2014-04-17T15:51:07 v2.6.0 User: admin Temperature: 37.0°C Pre 1 2 3 4



Literature

Haas T, Spielmann N, Mauch J, Madjdpour C, Speer O, Schmutz M, Weiss M. **Comparison of thromboelastometry (ROTEM) with standard plasmatic coagulation testing in paediatric surgery.** British Journal of Anaesthesia 2012 Jan; 108 (1): 36-4.

Theusinger OM, Nürnberg J, Asmis LM, Seifert B, Spahn DR. **Rotation thromboelastometry (ROTEM) stability and reproducibility over time.** Eur J Cardiothorac Surg. 2010 Mar; 37 (3): 677-83.


Haas T, Spielmann N, Mauch J, Speer O, Schmutz M, Weiss M. **Reproducibility of thromboelastometry (ROTEM): Point-of-care versus hospital laboratory performance.** Scand J Clin Lab Invest. 2012 Jul; 72(4):313-7.

Fast therapeutic decisions




- First results are available within 5 to 10 minutes
- 6 channels: differential diagnosis of coagulopathies
- Simplified interpretation of results via colour coded TEMograms/aggregation graphs and highlighted abnormal parameter results
- Easy therapeutic control using previous patient results as overlays
- Fast interpretation by overlay of standard curves over running results
- ROTEM® Data solutions enabled by HIS/LIS connection for comprehensive data transfer
- Mobile use via the ROTEM® trolley


ROTEM® analysis



CT	Clotting time	A6	Amplitude 6 min
CFT	Clot formation time	MS	Maximum slope
alpha	Alpha angle	AUC	Area under the curve
A10	Amplitude 10 min. after CT		
MCF	Maximum clot firmness		
LI30	Lysis index 30 min. after CT		
ML	Maximum Lysis		



TEMograms and aggregation curves of the 6 channels



Literature

Davenport R, Manson J, De 'ath H, Platten S, Coates A, Allard S, Pearse R, Pasi KJ, Maccallum P, Stanworth S, Brohi K. **Functional definition and characterization of acute traumatic coagulopathy.** Crit Care Med. 2011 Dec;36(12):2652-8.

Schöchli H, Maegele M, Solomon C, Görlinger K, Voelckel W. **Early and individualized goal-directed therapy for trauma-induced coagulopathy.** Scand J Trauma Resusc Emerg Med. 2012 Feb 24;20(1):15.


Görlinger K: **Gerinnungsmanagement bei Lebertransplantationen (Coagulation management during liver transplantation).** Hämostaseologie 2006; 26 (Suppl. 1): S64-S76.

Advanced diagnosis safety




- Instrument handling in a busy operating area enabled by the ball bearing stabilised technology of the ROTEM® thromboelastometry
- Differential diagnosis by the combination of up to 8 different assays
- The barcode scanner prevents the use of wrong or expired reagents
- Quality controls: ROTROL N (Level I); ROTROL P (Level II)
- Simple patient-ID search function for fast and safe real time data transmission

Barcode scanner



ROTEM® technology



Literature

Weber CF, Görlinger K, Meiniger D, Herrmann E, Bingold T, Moritz A, Cohn LH, Zacharowski K. **Point-of-care testing: A prospective, randomized clinical trial of efficacy in coagulopathic cardiac surgery patients.** Anesthesiology 2012 Sept; 117 (3): 531-547.


Larsen OH, Fenger-Eriksen C, Christiansen K, Ingerslev J, Sørensen B. **Diagnostic Performance and Therapeutic Consequence of Thromboelastometry Activated by Kaolin versus a Panel of Specific Reagents.** Anesthesiology 2011; 115:294-302

ROTEM® Data solutions



Data solutions - On-time information. Where it is needed.

All test results are available in numerical and graphical formats. Transfer to LIS/HIS systems or real time transmission to defined PCs is possible. Data protection via integrated password at the user management is ensured.



The complete ROTEM® system with ROTEM® *delta* and ROTEM® *platelet* provides an overview about the coagulation status within 10 minutes with information on:

- Requirement for factor, fibrinogen or platelet substitution
- Detection of platelet function and aggregation
- Hyperfibrinolysis
- Extent of dilutional coagulopathy
- Heparin and protamin dosage monitoring

India Contact

B-8, Sector-27
Noida - 201301
UP India

T: +91 (120) 4316140
F: +91 (120) 4275869

www.vijyoti.com | info@vijyoti.com
www.rotem.de | info@tem-international.de

